





Accessories for ULTRAMAX® DF and SF units

Create your own polishing Ferramentas by inserting Pedras and lapping Ferramentas into the appropriate steel clamp. Use the M6 (6mm thread) clamp Ferramenta Suportes in the UF-9700 ultrasonic handpiece supplied with ULTRAMAX DF and SF. For optimum machine performance, choose a medium clamp Ferramenta Suporte for Ferramentas up to 6.35mm in width and a Grande clamp Ferramenta Suporte for Ferramentas up to 15mm in width. Replacement screws

All drawing measurements are in millimeters



A. 6mm Thread Round Clamp Ferramenta Suporte

For round Ferramentas

For Dø (Diameter)	Item#
A. 2.0mm round Ferramenta	508-3508
2.35mm (3/32") round Ferramenta	508-3509
3.0mm round Ferramenta	508-3512*
3.3mm (1/8") round Ferramenta	508-3516
4.0mm round Ferramenta	508-3520

*Supplied with ULTRAMAX DF/SF

B. 6mm Thread Medium Clamp Ferramenta Suporte

For Plano Ferramentas up to 6.35mm (1/4") wide

For T (Thickness)	Item#
B. 0.5mm Plano Ferramenta	508-3495
0.8mm Plano Ferramenta	508-3500*
1.0mm Plano Ferramenta	508-3502*
1.5mm Plano Ferramenta	508-3504
2.0mm Plano Ferramenta	508-3506
3.0mm Plano Ferramenta	508-3510*
3.3mm (1/8") Plano Ferramenta	508-3514*

*Supplied with ULTRAMAX DF/SF

C. 6mm Thread Adjustable Clamp Ferramenta Suporte

For Plano Ferramentas 0.1mm to 8.0mm wide

For T (Thickness)	Item#
0.1mm to 1.5mm thick	508-3522

D. 6mm Thread Grande Clamp Ferramenta Suporte

For Plano Ferramentas 6.35mm to 15mm wide

For T (Thickness)	Item#
D. 1.0mm Plano Ferramenta	508-3503
1.5mm Plano Ferramenta	508-3505
2.0mm Plano Ferramenta	508-3507
3.0mm Plano Ferramenta	508-3511
3.3mm (1/8") Plano Ferramenta	508-3515
4.0mm Plano Ferramenta	508-3518

E, F ULTRAMAX® DF/SF Accessories and Replacements

Descrição	Item#
E. UF-9700 Handpiece only (6mm thread)	500-0205*
6/4mm Ferramenta Adapter (to hold M4 ultrasonic Ferramentas)	508-3342*
13.25	
Handpiece Cradle Rest	510-1130*
Standard Handpiece Protector Sleeve	500-0208*
F. Variable-Speed Foot Rheostat	500-0206